

CAPITAL CONNECT

INSTALL INTELEVIEWER

1. Open Internet –Google Chrome, Firefox, Explorer, Safari,
2. In the Address field type:
<https://connect.capitalradiology.com.au>
3. Select Download IntelViewer for the installer version that matches your desktop operating system (ie. 32 bit, 64 bit, Mac OS)
4. Double click the download installer file, and the installation process will begin
5. Accept the terms and press NEXT (next twoscreens)



6. Press NEXT (next 3 screens)



7. Click Done